



# 4K INTEGRATED MICROSCOPE

JBM-PL4



- ◆ HDMI HD output interface;
- ◆ One-click freeze/unfreeze function/on-screen comparison;
- ◆ Remote control focusing without shaking;
- ◆ Electronic zoom function;
- ◆ Supports manual/automatic exposure, and adjustment of exposure compensation value;
- ◆ Supports manual/automatic white balance, with adjustable red, green, and blue settings;
- ◆ Supports U-disk storage, with a maximum capacity of 64GB;
- ◆ Supports remote control operation, ensuring clearer images without shaking during photography;
- ◆ Supports 8 menu languages (default is English);
- ◆ Integrated design, featuring a built-in camera, light source, and microscope lens
- ◆ Supports saving shutdown settings parameters;



4K Resolution



HDMI Output



Remote Control



Freeze Function



Electronic Zoom



U-disk Storage



Built-in Camera



Built-in Light



## Specification

|                       |                                                                             |
|-----------------------|-----------------------------------------------------------------------------|
| Model                 | JBM-PL4                                                                     |
| sensor                | 1/2.8"                                                                      |
| Pixel size            | 1.45μm*1.45μm                                                               |
| Monitor               | 10.1-inch high-definition screen                                            |
| HDMI output           | 3840*2160 30fps 1920*1080 60fps                                             |
| USB output            | 3840*2160/2560*1440<br>1920*1080/1280*720<br>640*480/424*240 320*240 (MJPG) |
| Packing size          | 315*145*275mm                                                               |
| Weight (gross weight) | 2.05 KG                                                                     |

# Products informations





J-9/36

RAJOURIGARDEN

NEW DELHI 110027

[www.jbmicroscope.com](http://www.jbmicroscope.com)

[info@jbinstruments.com](mailto:info@jbinstruments.com)

+91 9354741245