

500x Upright Microscope

JBM-UM50



[illegible]

Product Features

Brand Name	JB MICROSCOPE
MODEL	JBM-UM50
Stage	140mmX120mm double-layer movable platform
Stage moving range	80mmX50mm
Focusing	Coarse and fine focusing, using gear rack transmission mechanism, fine-motion value 0.002mm
Condenser	NA1.25 swing Abbe condenser with variable aperture
Lighting source	Upper lighting LED lamp SW, DC12V2A, brightness adjustable
Down lighting	LED lightSW, DC12V2A, adjustable brightness
Light collecting mirror	high brightness fixed lighting with adjustable light bar

Product Technical Parameter

Objective lens (infinity plan achromatic Stress-free metallographic objective lens)

Magnification	Numerical Aperture	Working Distance (mm)
5x 10x 20x 50	0.12 0.25 0.4 0.55	23.6,17.7,10.4,7

Total Magnification after the combination of eyepiece and objective lens:

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++