

ESD Stereo Zoom Microscope JBM-EPL



Protect your sensitive components from electrostatic discharges with JBM ESD Safe Trinocular Microscope. Provides crisp clear 3D images, long working distance and a wide field of view. 0.67x to 45x continuous zoom with 10x eyepieces. View magnified images through the microscope eyepieces or on a monitor.





Technical Specifications

Brand	JB MICROSCOPE
MODEL NO	JBM-EPL
Magnification	45x
Eyepiece	WF10X/22MM high eyepoint, wide-field; Convenient for observers wearing glasses.
Objective	0.67X-4.5X Stereo Zoom
Zoom Ration	6:3:1
Working Distance	150mm
Stand Type	Boom Stand
Vertical Post Height	380mm
Maximum Vertical Post Extended Length	254mm
Base Type / Shape	Heavy Duty Base, Rectangle
Base Dimension	285x260x15mm
Surface Treatment	Electroplating Black / Silver



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com | info@jbinstruments.com | +91 9354741245