

ESD Stereo Zoom Microscope

JBM-PL45



[illegible]

[illegible][illegible][illegible][illegible][illegible][illegible][illegible][illegible]

Technical Parameters

Brand	JB Microscope
Viewing Head	45° inclined binocular head with fixed eyepiece tube, inter-pupillary range: 52mm – 76mm, with click stop and interlock mechanism
	45° inclined trinocular head with fixed eyepiece tube, inter-pupillary range: 52mm – 76mm, with click stop and interlock mechanism
Eyepiece	High eye-point wide field plan eyepiece PL10X/22mm, with adjustable diopter
Main Objective	Zoom objective 0.67X – 4.5X,
working distance	110mm
Auxiliary Objective	0.5X/177mm, 0.7X/120mm, 1.5X/47mm
Stand	Arm stand with fixed focus holder
	Pillar stand, can be equipped with different focus holders
Base	Plan base
	Base with illumination
Camera adapter	0.35X/0.5X/0.65X/1X C-mount adapted
Others	Precision measuring stage, gem clip

OPTIONAL

ESD SAFETY

JBМ-PL-45 Series microscopes feature a specialized ESD Safe coating, protecting sensitive components from electrostatic discharges

OPTIONAL ACCESSORIES

LCD MicroscopeDigital Camera

JBM-L13



13-INCH

LCD is a multifunctional all-in-one camera that integrates a 13.3-inch display screen and a 4K image acquisition and processing system for capturing digital images of stereo microscopes, biological microscopes, or online interactive teaching.

DIGITAL MICROSCOPE CAMERA

JBM-230-4K



HIGH PRECISION INDUSTRIAL SMART CAMERA
JBM SERIES SDK 60FPS INDUSTRIALFULLY
-INTEGRATED FLEXIBILITY OCR.

FHD SCREEN



24-inch

HD screen for crystal-clear viewing
Full HD Display Brings Microscopic
Images To LifeWith Exceptional
Color And Contrast.

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com | info@jbinstruments.com | +91 9354741245

+++++