

# ESD Stereo Zoom Microscope JBM-EPLT





# ESD Safe Trinocular Microscope

Protect your sensitive components from electrostatic discharges with JBM ESD Safe Trinocular Microscope.

Provides crisp clear 3D images, long working distance and a wide field of view. 0.67x to 45x continuous zoom with 10x eyepieces. View magnified images through the microscope eyepieces or on a monitor.



## Crisp Clear 3D Images

High resolution, sharp and detailed visuals.



## Long Working Distance

More space for comfortable and efficient work.



## Wide Field of View

See more, work faster.



## 0.67x to 45x Continuous Zoom

With 10x eyepieces.



## Dual Viewing

View magnified images through eyepieces or on a monitor.



SPECIFICATION

|                                       |                                                                                 |
|---------------------------------------|---------------------------------------------------------------------------------|
| Brand Name                            | JB MICROSCOPE                                                                   |
| MODEL                                 | JBM-EPLT                                                                        |
| Magnification                         | 45x                                                                             |
| Eyepiece                              | WF10X/22MM high eyepoint, wide-field; Convenient for observers wearing glasses. |
| Objective                             | 0.67X-4.5X Stereo Zoom;                                                         |
| Zoom Ration                           | 6:3:1                                                                           |
| Working Distance                      | 150mm                                                                           |
| Stand Type                            | Boom Stand                                                                      |
| Vertical Post Height                  | 380mm                                                                           |
| Maximum Vertical Post Extended Length | 254mm                                                                           |
| Base Type / Shape                     | Heavy Duty Base, Rectangle                                                      |
| Base Dimension                        | 285x260x15mm                                                                    |
| Surface Treatment                     | Electroplating Black / Silver                                                   |

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

[www.jbmicroscope.com](http://www.jbmicroscope.com)

[info@jbinstruments.com](mailto:info@jbinstruments.com)

+91 9354741245

+++++