

Horizontal Digital Inspection Microscope JBM-HZ45i




Features:

- Zoom magnification power offering sharp images
- Metal framework with solid base plate
- 0.7X-5.0X zoom objective
- C-mount lens
- Long working distance
- Sharp images over a super wide field view



Product Specifications

Brand	JB Microscope
Microscope	Video Microscope
Drawtube	Monocular
Material	Aluminium Alloy
Application	PCB Inspection
Stage	Plain Stage
Easy Large	Sized Focusing Knobs on Both Sides of Focusing Holder.Steel Gears on Focusing Holder.
Standard Size Focusing Holder	Will Fit Scope Heads with Diameter of 50mm
Base size	380mm X 250mm
Frame diameter	50mm
Load table size	180mm X 120mm
Load table adjustable stroke	X: 70mm, Y: 40mm



+++++

+++++

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

WWW.JBM ICROSCOPE.COM

INFO@JBINSTRUMENTS.COM

+91 9354741245

+++++