






LCD MICROSCOPE DIGITAL CAMERA JBM-L10



- New generation 10-inch large display with a built-in camera, built-in WiFi transmitter.
- The mobile phone can be connected to WiFi, synchronous observation.
- 10-inch high-definition large screen, the image is clearer.
- You can rotate it by 90 degrees.
- The ultra-long working distance is conducive to spot welding by mobile phone maintenance, circuit board inspection, etc.
- Remote button remote control / External fill light / Free reflector.
- Continuous zoom trinocular simultaneous observation, clear image, easy operation.
- Large wide-angle lens / continuous zoom / adjustable LED / diopter adjustment.



Specification

Brand	JB MICROSCOPE
Model	JBM-L10
Screen Size	10 Inch HD LCD Display
Resolution	1080P
Brightness	Adjustable
Interface	Micro / USB 2.0, HDMI
Power Supply	3000mAh
Memory Card	16 GB TF Card Include (Max support 128G card)
Output Port	HDMI
Viewing Angle	IPS Full View

+++++

+++++

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++