



MICROSCOPE AUTOFOCUS CAMERA

JBM-AF8



Autofocus

Fast & Accurate
Focusing



High Resolution

Crystal Clear
Images



Multiple Output

HDMI&USB
Interface



Easy to Use

Plug & Play
Setup



PRECISE AUTOFOCUS

Consistent & Reliable Results



ADVANCED SENSOR

Superior Image Quality




BUILT FOR RELIABILITY

Durable. Stable. Professional.

JBM-AF8

JBM-AF8 Plus can measure straight lines, angles, concentric circles, parallel lines, point-to-line distance, point-to-line-to-concentric circles, diameter, circumference and area, length and width of rectangle. It can measure the polygon and calculate its area. All the corresponding measurement report can be output to U disk.


Find of View



0.7x Magnification



1x Magnification



2x Magnification




Technical Specification

Brand	JB Microscope
Image Sensor	Color CMOS
Pixel	8 MI pixels
Pixel size	1/2.8"
Menu	All-digital UI design
Method of Operation	Mouse
Power DC	DC12V
output method	HDMI
Display Frame rate	4K@60fps/4K@50fps
Scanning method	Line by line scanning
White Balance	Manual / Auto/ One click
Exposure	Auto/Manual
Shutter Speed	1/50s(1/60s)~1/10000s
Angle of view	D:78.6°(W)-6.16°(T) H:69.2°(W)-5.22°(T) V:45.4°(W)-3.26°(T)
Working distance	10cm~∞
Optical magnification	1~28
Digital magnification	1~3
Digital zoom	Support
Storage	Support USB2.0/USB3.0 external storage Support camera internal storage


SOFTWARE INTERFACE DISPLAY

Measuring tools for point, lines ,polylines, rectangles ,custom polygons, parallel lines, angles, concentric circles, center distance, two -point fixed circles,three-point fixed circles,center lines, crosshairs,arcs,wire spacing and more

Introduction to the main interface of the monitor



Introduction to the main interface of the laptop



Operating system: Windows XP 32bit ;Windows 7/8/10 32/64bit

OPTIONAL ACCESSORIES

4KFHD SCREEN



27-inch

4K Ultra FHD screen for crystal-clear viewing. Full HD Display Brings Microscopic Images To Life With Exceptional Color And Contrast.

2KFHD SCREEN



24-inch

HD screen for crystal-clear viewing. Full HD Display Brings Microscopic Images To Life With Exceptional Color And Contrast.

SCREEN



17-inch

Screen for crystal-clear viewing. Display Brings Microscopic Images To Life With Exceptional Color And Contrast.



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245