

Stereo Zoom Microscope JBM-765T



Technical Specifications

Brand Name	JB MICROSCOPE
MODEL	JBM-765T
Trinocular Head	45° inclined 360°rotate
Interpupillary Distance Adjustment	54-76mm
Diopter Adjustment	±5dp
Eyepiece	WF10X/22mm
Objective	0.65X-6.5X
Magnification Ratio	1:10
Magnification	0.65X-65X
Working Distance	110mm
C-mount Adapter	1x
Vertical Post Height	240mm
Vertical Post Diameter	φ32mm
Microscope Body Holder	φ76mm
Focus Travel	50mm
Stage Plate	φ95mm black white contrast plate
Base	274x207x23mm
Material	Metal

[illegible]

+++++



J-9/36

RAJOURIGARDEN

NEW DELHI 110027

www.jbmicroscope.com

info@jbinstruments.com

+91 9354741245

+++++